

NEW OKOLAB® INCUBATING CHAMBERS TO ENABLE LIVE CELL IMAGING
ON THE PRIOR® QUEENSGATE SP200, SP400, SP450 AND SP600
NANOPOSITIONING PIEZO SAMPLE SCANNERS



On stage gas chamber for
Okolab Cage Incubators



ORDER CODE: **H201-PRIOR-SP200-400-600**

On stage heated chamber for
Okolab UNO and Bold Line
Top Stage Incubators



ORDER CODE: **H301-PRIOR-SP200-400-600**


FEATURES:

- Easy coupling to the Prior piezo
- Compatible with both R and XR Prior motorized stages as well as with Nikon XY motorized stages underneath the piezo
- Great variety of multi-accommodation interchangeable magnetic sample holders
- Compatible with perfusion applications
- Compatible with sample feedback thermal sensor

CAN ACCEPT ANY TYPE OF SAMPLE VESSEL, PLEASE FIND BELOW COMPATIBLE FORMATS. MULTI-ACCOMMODATION HOLDERS' MODELS ARE AVAILABLE.




35mm petri dish




1"x3" slide




60mm petri dish



1"x2" Lab-Tek chambered cover glass type I and type II



1"x3" chamber slide



6-12-24-48-96 well plates

MICROSCOPE INCUBATOR MODELS

CAGE INCUBATOR



- Microscope enclosure
- Obscuring panels (*optional*)
- Digital temperature control
- Filtering box (*optional*)
- Digital gas control
- Humidity module
- Stage top gas chamber
- Operated via touch screen interface with on-board memory
- Compatible with Smart Box for remote control

STAGE TOP INCUBATOR



- Digital temperature control
- Digital Gas controller (Bold Line) or Manual Gas Mixer (UNO)
- Digital humidity module
- Operated via touch screen interface with on-board memory
- Compatible with Smart Box (Bold Line) for remote control



okolab

www.oko-lab.com

Okolab Srl

Via A. Olivetti, 1
80078 Pozzuoli, NA, Italy

Luca Lanzaro, Ph.D

+39 081 8062624

lanzaro@oko-lab.com

Okolab USA Inc.

233 Merchant Street, Suite 500
Ambridge, PA 15003, USA

Silvia Foppiano, Ph.D

+1 (650) 228 6154

foppiano@oko-lab.com

Lara Petrak

+1 (650) 410 0756

petrak@oko-lab.com

Okolab Shanghai Co., Ltd.

2F No.388 Madang Rd.,
Premise Suite #37, Huangpu Dist.,
Shanghai 200020, China

Pony Tong

+86 18500508458

tong@oko-lab.cn